




Y-Griff ohne Sperre
Axial handle without
ratchet



Y-Griff mit Sperre
Axial handle with
ratchet

Ø 2,7 mm



27970-01

Universalfaßzange, 135 mm AL

Universal grasping forceps, 135 mm wl

27970-21 = as above but 250 mm wl

27970-71 = as above but 350 mm wl




27971-01

Atraumatische Faßzange, 135 mm AL

Atraumatic grasping forceps, 135 mm wl

27971-21 = as above but 250 mm wl

27971-71 = as above but 350 mm wl



27971-07

Atraumatische Allisklemme, 135 mm AL

Atraumatic Allis-clamp, 135 mm wl

27971-27 = as above but 250 mm wl

27971-77 = as above but 350 mm wl



27971-08

Atrauma Organfaßzange, 135 mm AL

Atraumatic Organ grasping forceps, 135 mm wl

27971-28 = as above but 250 mm wl

27971-78 = as above but 350 mm wl



27971-10

Gefensterte Faßzange, 135 mm AL

Fenestrated grasping forceps, 135 mm wl

27971-30 = as above but 250 mm wl

27971-80 = as above but 350 mm wl



27971-25

Atraumatische Pinzette, 135 mm AL

Atraumatic dressing forceps, 135 mm wl

27971-75 = as above but 250 mm wl

27971-85 = as above but 350 mm wl



27973-01

Dissektionsklemme, 135 mm AL

Dissection clamp, 135 mm wl

27973-21 = as above but 250 mm wl

27973-71 = as above but 350 mm wl



27973-13

Atrauma Faßzange X-cross, 135 mm AL

Atraumatic grasping forceps x-cross, 135 mm

27973-53 = as above but 250 mm wl

27973-83 = as above but 350 mm wl



27974-01

Hakenschere, 135 mm AL

Hook scissors, 135 mm wl

27974-21 = as above but 250 mm wl

27974-71 = as above but 350 mm wl



27975-01



Schere gerade, 135 mm AL

Scissors curved, 135 mm wl

27975-21 = as above but 250 mm wl

27975-71 = as above but 350 mm wl



27975-02



Schere gebogen, 135 mm AL

Scissors curved, 135 mm wl

27975-22 = as above but 250 mm wl

27975-72 = as above but 350 mm wl



27976-01



Biopsiezange, 135 mm AL

Biopsy forceps, 135 mm wl

27976-21 = as above but 250 mm wl

27976-71 = as above but 350 mm wl



27968-05



Top dur Nadelhalter gerade, 135 mm AL

Top dur Needle holder straight, 135 mm wl

27968-25 = as above but 250 mm wl

27968-75 = as above but 350 mm wl



27968-06



Top dur Nadelhalter gebogen, 135 mm AL

Top dur Needle holder curved, 135 mm wl

27968-26 = as above but 250 mm wl

27968-76 = as above but 350 mm wl



27968-35



DIMICRON Nadelhalter konvex/konkav,
135 mm AL

DIMICRON Needle holder convex/concave,
135 mm wl

27968-85 = as above but 250 mm wl

27968-65 = as above but 350 mm wl



27968-09




Top dur Gegennadelhalter gerade, ohne
Sperre, 135 mm AL

Top dur Revers needle holder straight, without
ratchet, 135 mm wl

27968-29 = as above but 250 mm wl

27968-79 = as above but 350 mm wl



27982-51



Retraktor, 250 mm AL

Retractor, 250 mm wl



**Pistolengriff und Pinzettengriff auf
Anfrage lieferbar**



**Pistol and pencil grip available on
request**





DUFNER Universal Retraktorsystem

- mit 16 möglichen Haltepositionen der Fixierungsvalven
- verschiedene Valvenvarianten: rund, Standard, vergrößert
- alle Valven mit leicht aufgebogem Ende für verbessertes Handling
- Ovalform für maximale Einsatzmöglichkeiten
- DIRONDAL Design = gerundete Formen aller Komponenten bis ins Detail ermöglichen eine maximal atraumatische Anwendung
- 3 Ringgrößen: Neonatal, Kleinkind und Jugendalter
- Verwendung von streng selektiertem deutschen Instrumentenstahl für höchstmögliche Stabilität
- feinste, präzise Verarbeitung für optimales Handling und besonders leichtes Gesamtgewicht

DUFNER Universal retractor system

- with 16 possible holding positions of the fixation valves
- different variations of valves: round, standard, enlarged
- all valves with slightly curved up ends for an improved handling
- oval form for maximal operation potential
- DIRONDAL design = round shaped forms of all components down to the last detail enables a maximal atraumatic usage
- Three ring sizes: neonatal, small child and adolescence
- Using of strict selected German instrument stainless steel for highest possible stability
- Finest precise finishing for optimal handling and special light total weight



12881-21

Ovalrahmen allein X130 x Y100 mm für
DUFNER Universal Retraktorsystem

Ovalform Frame only X130 x Y100 mm
for DUFNER Universal Retractorsystem

12881-22

Ovalrahmen allein X100 x Y75 mm für
DUFNER Universal Retraktorsystem

Ovalform Frame only X100 x Y75 mm
for DUFNER Universal Retractorsystem

12881-23

Ovalrahmen allein X75 x Y57 mm für
DUFNER Universal Retraktorsystem

Ovalform Frame only X75 x Y57 mm
for DUFNER Universal Retractorsystem

12881-31

Runde Valve 16 x 12 mm für Ovalrahmen
12881-21/-22/-23

Round valve 16 x 12 mm for ovalform frame
12881-21/-22/-23

12881-32

Standardvalve 25 x 8 mm für Ovalrahmen
12881-21/-22/-23

Standard valve 25 x 8 mm for ovalform
frame 12881-21/-22/-23

12881-32-000

Standardvalve 60 x 25 mm für Ovalrahmen
12881-21/-22/-23

Standard valve 60 x 25 mm for ovalform
frame 12881-21/-22/-23

12881-32-001

Standardvalve 60 x 16 mm für Ovalrahmen
12881-21/-22/-23

Standard valve 60 x 16 mm for ovalform
frame 12881-21/-22/-23

12881-32-005

Standardvalve 60 x 25 mm, C-Valvenform
für Ovalrahmen 12881-21/-22/-23

Standard valve 60 x 25 mm, c-form for
ovalform frame 12881-21/-22/-23

12881-32-006

Standardvalve 60 x 16 mm, C-Valvenform
für Ovalrahmen 12881-21/-22/-23

Standard valve 60 x 16 mm, c-form for
ovalform frame 12881-21/-22/-23

12881-33

Standard Valve 30 x 16 mm für Ovalrahmen
12881-21/-22/-23

Standard valve 30 x 16 mm for ovalform
frame 12881-21/-22/-23