



### Laparoscopic Needle Holders

<b>PRODUCT</b>	Laparoscopic Needle Holders
<b>ARTG NO:</b>	270744

#### PRODUCT DESCRIPTION

Comprehensive range of Laparoscopic Needle Holders. These high-quality instruments are German designed and manufactured.

#### FEATURES:

- Choice of handle design: with aluminium anodised axial handle – straight or curved, with stainless steel axial handle
- Ergonomic design
- Straight, curved left, curved right and self-righting jaw options
- Tungsten Carbide inserts optional
- 5mm diameter, 33cm and 45cm lengths
- Affordable
- Re-usable
- Easy cleaning and processing

#### PRODUCT IMAGES

