

Liquiband Fix8 Open

PRODUCT	Open Hernia Mesh Fixation Device
PRODUCT CODES	ML-FX002
ARTG No	314884
BILLING CODE	WL001

PRODUCT DESCRIPTION

The LiquiBand FIX8® Open Hernia Mesh Fixation device is designed for the application of cyanoacrylate adhesive to an implanted hernia repair mesh, in order to fix the underlying tissue.

Designed for atraumatic hernia mesh fixation in open surgery.

FEATURES:

- Strong and secure mesh fixation
- Precise and controlled application
- Fixation at multiple angles
- New atraumatic tip design
- Enhanced anchor delivery
- Peritoneum closure
- Indicated for ventral and incisional hernia
- Sterile, single-use

PRODUCT IMAGE



PACKAGING SPECIFICATIONS

• Individually Packed	• Box of 1 UOM	• Sterile Single Use
-----------------------	----------------	----------------------