

Panther Manual Endo Stapler	
PRODUCT	PANTHER ENDO Linear Cutter Stapler, with non-stop Articulation and grasping mode
PRODUCT CODES	CEAC-30, CEAC-45, CEAC-60
ARTG No	389982

PRODUCT DESCRIPTION

The Endo Linear cutter stapler, single-use is intended for transection, resection, and/or creation of anastomoses. The instruments have application in multiple open or minimally invasive general, gynaecological, colorectal, bariatric and thoracic surgical procedures.

FEATURES:

- Innovative 45 degrees non-stop articulation system
- Universal handle ready to use with most extended portfolio reloads
- Ergonomic universal handle design with 360-degree rotation
- Grasping tissue mode
- Can be used with P3H and P2G Reload technology
- Easy to assemble and use
- Sterile, single-use

PRODUCT IMAGE



PACKAGING SPECIFICATIONS

• Individually Packed	• Box of 1 UOM	• Sterile Single Use
-----------------------	----------------	----------------------