

## Single-Use BiPolar Forceps

We offer a complete range of disposable bipolar forceps, all with integrated cables and options of standard or non-stick tips

### Features

- Lightweight and quality finish single-use BiPolar Forceps
- The comprehensive range consists of all styles and sizes including McPherson, Jeweller, Adsons, Straight and Bayonet. 85mm to 250mm lengths - 0.5mm, 0.7mm and 1mm straight and angled tips
- Packaged as 1pc/pouch, 10 pouches/box
- Sterile, single use
- Blue – Standard tips, Pink – Non-stick tips

**ARTG No:** 362852

