



Laparoscopic Maryland Bipolar Forceps

PRODUCT	Laparoscopic Maryland Bipolar Forceps
PRODUCT CODES	109Y.202
ARTG No	365071

PRODUCT DESCRIPTION

The Laparoscopic Maryland Bipolar Forcep is designed to dissect and coagulate selected tissue in laparoscopic surgery. It is to be connected through a suitable cable with the bipolar output of an electro-surgical unit/generator.

FEATURES:

- For secure grasping and precise dissection
- High level of tactility
- Excellent coagulation results
- Safe bipolar technology
- Complete with integrated fixed-width pin cable
- Double action
- 5mm diameter x 33cm length
- Sterile, single-use
- Box 5

PRODUCT IMAGE



PACKAGING SPECIFICATIONS

• Individually Packed	• Box of 5 UOM	• Sterile Single Use
-----------------------	----------------	----------------------