



Australian Designed Safer Delivery System



VARIVAC™

THUMB GRIP TROCAR

- The Patented "Thumb Grip Trocar" gives the Surgeon more control which helps prevent needlestick injuries during insertion.

SUPERIOR DRAINAGE

- Four flow channels for greater tissue contact. No holes means less tissue ingrowth and less risk of occlusion.

PATIENT COMFORT

- Silicone tubing allows for more flexibility.
- The design means less tissue trauma on removal of drain.
- Latex Free

FEATURES





- Radio opaque.
- Available in round, flat and multiple channel.
- With or without trocar.
- All drains are double wrapped sterile, with a connector that will suit most types of wound drainage systems.
- Compatible with high and low suction.
- Visible skin mark.






Australian Designed, Latex Free Wound Drains



Ordering Information **Without** Trocar

	Code	Description	Size	Case Qty
Round Silicone				
	VC500-32	Round Silicone Drain (with side holes)	10Fr	10
	VC500-33	Round Silicone Drain (with side holes)	16Fr	10
	VC500-34	Round Silicone Drain (with side holes)	20Fr	10
Flat Silicone - Single Channel				
	VC500-37	Flat Silicone Drain (with side holes) 3/4 length holes in drain section	7mm	10
	VC500-38	Flat Silicone Drain (with side holes) 3/4 length holes in drain section	10mm	10
VariBlake™ Silicone - Multiple Channel				
	VC500-41	VariBlake™ Silicone Drain without trocar (390mm channel section & 800mm tube without perforation)	10Fr	10
	VC500-42	VariBlake™ Silicone Drain without trocar (390mm channel section & 800mm tube without perforation)	15Fr	10
	VC500-43	VariBlake™ Silicone Drain without trocar (390mm channel section & 800mm tube without perforation)	19Fr	10
	VC500-46	VariBlake Silicone Drain without trocar (350mm channel section)	24Fr	10
VariBlake™ Silicone - Multiple Channel Flat Drain				
	VC500-44	Flat VariBlake™ Silicone Drain without trocar (390mm perforation & 800mm drain without perforation)	8mm	10
	VC500-45	Flat VariBlake™ Silicone Drain without trocar (390mm perforation & 800mm drain without perforation)	11mm	10

Ordering Information **With** Trocar

	Code	Description	Size	Case Qty
VariBlake™ Silicone - Multiple Channel				
	VC500-51	VariBlake™ Silicone Drain with trocar (390mm channel section & 800mm tube without perforation)	10Fr	10
	VC500-52	VariBlake™ Silicone Drain with trocar (390mm channel section & 800mm tube without perforation)	15Fr	10
	VC500-53	VariBlake™ Silicone Drain with trocar (390mm channel section & 800mm tube without perforation)	19Fr	10
	VC500-53ST	VariBlake™ Silicone 19Fr Drain with 16Fr trocar (390mm channel section & 800mm tube without perforation)	19Fr	10
VariBlake™ Silicone - Multiple Channel Flat Drain				
	VC500-54	Flat VariBlake™ Silicone Drain with trocar (390mm perforation & 800mm drain without perforation)	8mm	10
	VC500-55	Flat VariBlake™ Silicone Drain with trocar (390mm perforation & 800mm drain without perforation)	11mm	10
T Drain				
	VC500-71	Silicone T Drain 10Fr drain only	10Fr	10
	VC500-72	Silicone T Drain 14Fr drain only	14Fr	10

™ * Trademark of International Medical Research.



Ascend Surgical and Medical
Unit 9/3 Leo Alley Road
Noosaville, Queensland 4566
Contact: 1300 784 893

Contact the team at Ascend Surgical and Medical for support
with the VariVac™ range of products.

www.ascendsurgical.com.au