

## Panther PSMART PLUS Endo Stapler

PRODUCT	PANTHER PSMART PLUS Fully Powered Endoscopic cutting stapler, Single-use
PRODUCT CODES	PEPS, PEPM, PEPL
ARTG No	386366

### PRODUCT DESCRIPTION

The PSMART PLUS Fully Powered endoscopic single use cutting stapler is intended for transection, resection, and/or creation of anastomoses. The stapler has application in multiple open or minimally invasive general, colorectal, bariatric and thoracic surgical procedures.

### FEATURES

- Fully Powered Rotation and Articulation
- “One Touch” System for One Handed Use
- SMART articulation control
- SMART reload identification
- SMART tissue recognition
- SMART firing control
- 15 Firings
- Universal fitting to Panther reloads, can be used with P3H and P2G Reload technology
- Sterile, single-use

### PRODUCT IMAGES



### PACKAGING SPECIFICATIONS

• Individually Packed	• Box of 1 UOM	• Sterile Single
-----------------------	----------------	------------------