

## Endoscopic Clip Appliers

|               |                                                                                                                                                                                                                                                         |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| PRODUCT       | InnoClip Endoscopic Clip Appliers<br>InnoClip Disposable Clip Applier, 330mm*5mm, Titanium, ML, 16 clips<br>InnoClip Disposable Clip Applier, 350mm*10mm, Titanium, L, 20 clips<br>InnoClip Disposable Clip Applier, 300mm*10mm, Titanium, ML, 20 clips |
| PRODUCT CODES | DC5T1-B6, TDC-10-L, TDC-10-ML                                                                                                                                                                                                                           |
| ARTG No       | 461225                                                                                                                                                                                                                                                  |

### PRODUCT DESCRIPTION

Ascend Surgical and Medical's Endoscopic Clip Appliers for use in endoscopic/laparoscopic procedures, to ligate vessels and tubular structures.

### FEATURES

- Multi fire with titanium clips for continuous ligation
- Audible and tactile feedback
- Anti-slip trigger
- Clip indicator window
- Easy one handed 360° rotation
- Inner cross textured clip for firm grip on vessels
- Advanced firing mechanism
- Affordable
- Individually packaged, Box 5
- Sterile, Single-Use
- Sterilisation - Ethylene Oxide

### PRODUCT IMAGES

