

Panther CADL 8mm Reload Unit

PRODUCT	Endo Articulating reload unit Curve Tip 8 mm diameter Articulated Reload - 45mm size with Fixed Anvil. 2 rows of staples - Diameter 8 mm - Staple height 2.5 mm - White
PRODUCT CODES	CADL-45T
ARTG No	386361

PRODUCT DESCRIPTION

The New CADL 8mm Reload Unit has a shaft length of 240mm and a total length of 330mm, is a 45mm reload with 4 rows of staples and with a 45° Angle. The CADL 8mm Reload Unit is for the following applications:

- Thoracic surgery for vascular anastomosis
- Vessels with difficult access - (Segmentectomy VATS -Uniportal VATS)
- Live donor kidney transplantation

FEATURES

Offering True 8mm Access with a narrow 6 mm anvil and 8 mm shaft, the CADL small diameter reload:

- Improves ability to gain access
- Aids in dissection around delicate tissue structures
- Facilitates vessel placement
- Compatible with all Panther Smart Endo series and Manual series endo staplers
- New knife for each reload unit

PRODUCT IMAGES



PACKAGING SPECIFICATIONS

- | | | |
|-----------------------|----------------|----------------------|
| • Individually Packed | • Box of 6 UOM | • Sterile Single Use |
|-----------------------|----------------|----------------------|