

Veress Needles

PRODUCT	Veress Needle 2.1mm*120mm and 150mm lengths
PRODUCT CODES	105Y.202, 105Y.205
ARTG No	469429

PRODUCT DESCRIPTION

Ascend Surgical and Medical Veress Needles for creating pneumoperitoneum.

FEATURES

- Sharp tip for easy insertion
- Spring loaded stylet provides audible 'click' when passing through tissue layers
- Red visual needle position indicator
- High-flow stopcock
- Affordable
- Individually packaged
- Sterile, Single-Use
- Box 10
- Sterilisation - Ethylene Oxide

PRODUCT IMAGES

