

Smart Line Single-Use Bipolar Forceps

PRODUCT	Smart Line Single-Use BiPolar Forceps, Single-use
PRODUCT CODES	STANDARD TIP / SLNA102R-05, SLNA102R-10, SLNA102R-10SR, SLNA131R-05, SLNA101R-10, SLNA108A-10, SLNA108A-05, SLNA118R-10, SLNA108R-15, SLNA125R-05, SLNA125R-10, SLD1303R-05, SLD1303R-10
	NON STICK TIP / SLNS102R-05, SLNS102R-10, SLNS102A-10, SLNS131R-05, SLNS108A-05, SLNS108A-10, SLNS118R-10, SLNS125R-05, SLNS125R-10, SLD1303RNS-05, SLD1303RNS-10
ARTG No	362852

PRODUCT DESCRIPTION

We offer a complete range of disposable bipolar forceps, all with integrated cables and options of standard or non-stick tips.

FEATURES

- Lightweight and quality finish single-use BiPolar Forceps
- The comprehensive range consists of all styles and sizes including McPherson, Jeweller, Adsons, Straight and Bayonet. 85mm to 250mm lengths, 6 and 8mm straight and angled tips
- Packaged as 1pc/pouch, 10 pouches/box
- Sterile, single use
- Blue – Standard tips, Orange – Non-stick tips

PRODUCT IMAGES



PACKAGING SPECIFICATIONS

- | | | |
|-----------------------|----------------|------------------|
| • Individually Packed | • Box of 1 UOM | • Sterile Single |
|-----------------------|----------------|------------------|